hydro surgery

Ausgabe 7

Das neue Cellgraft®-System

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280 patients were treated from October 1999 until March 2003."

Results: "The targeted fragmented removal of the fat tissue allows the achievement of the desired result in a controlled and safe way. Surgery times can be reduced by more than 40 %. The technique of water-jet liposcution is a safe, gentle and targeted method to remove subcutaneous fat build-ups. It offers a very good way for molding the tissue during the operation. Drug-related side effects are not to be expected. The method is simple, easily explained to the patient and quickly learned by the surgeon. The postoperative leakage of fluid from the incisions with the tumescent method is largely reduced if not even stopped with the new water-jet method."

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Araco, A., Gravante, M.D., Araco, F., Delogu, D., Cervelli, V.: Comparison of Power Water-Assisted and Traditional Liposuction: A Prospective Randomized Trial of Postoperative Pain. Aesth. Plast. Surg. 31: 259 – 265: 2007.

A prospective randomized trial compared power water-assisted liposuction with the traditional tumescent technique. Patients were randomly assigned to water-assisted or traditional liposuction.

Results: "From September 2005 to December 2005, 60 patients were recruited and analyzed. For the study, 28 patients were randomized to traditional liposuction and 32 to power water-assisted liposuction. A significant difference in postoperative pain was observed (p < 0.05). After 4 days, 87 % of the patients who underwent power water-assisted

The study findings demonstrate that water-assisted liposuction is an almost painless procedure that produces less tissue trauma than traditional liposuction."

liposuction were completely pain free, as compa-

red with 3.6 % of those treated with traditional

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Man, D.; Meyer, H.: Water Jet-Assisted Lipoplasty. Aesthetic Surgery Journal; May/June 2007, 342 – 346.

In total, the authors "have carried out more than 800 treatments since 2001".

Results: "Compared with the quantity of tumescent solution used in conventional manual lipoplasty (100 %), with WAL only an average of 20 % to 30 % was used in preinfiltration. The length of time tumescent solution remains in the tissue, as well as the resulting absorption times, are all considerably lower for the recommended infiltration solutions compared with all other tumescence-based lipoplasty techniques." "The average patient satisfaction rate for both the procedure and the final outcome was higher than 94 %."

Conclusions: "On the basis of our clinical experience with body-jet WAL under local anesthesia, we have found that:

- 1. WAL presents a new and safe method that is suitable for all types of lipoplasty.
- 2. In almost all cases, WAL facilitates use of preinfiltration of modified tumescent solution to create analgesia that is suitable for the performance of painless or near painless lipoplasty.
- 3. There is significantly reduced painrelated impairment during and after the procedure compared with standard tumescent technique. Patients recover quickly and return to normal daily activities rapidly.
- 4. Considerably less intraoperative swelling allows the surgeon to realize the target result with greater precision
- 5. Fine shaping of small deposits with precision is also possible."

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Stutz, J.J. D. Krahl: Water-Jet Assisted Liposuction for Patients with Lipoedema: Histologic and Immunohistologic Analysis of the Aspirates of 30 Lipoedema Patients. Aesthetic Plastic Surgery (2009)33: 153-162.

"WAL liposuction treatment can bring longterm improvement if the operative technique focuses on lymph vessel preservation. Immunohistologic analyses show minimal evidence of lymph vessel structures in lipoaspirates. The histologic analysis of the aspirates documents a relatively specific removal ("apheresis") of primarily intact lipocytes with low vascular amount."

"After Water-Assisted Liposuction with the bodyjet®, the lipocytes in the aspirate were shown to be predominantly intact: "In 28 of the 30 investigated lipoaspirates (patients), the lipocytes were found to be predominantly (>70 %) intact."

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G.H. Sasaki: Water-Assisted Liposuction for Body Contouring and Lipoharvesting - Safety and Efficacy in 41 Consecutive Patients. Aesthetic Surgery Journal 2011; 31: 76.

"Fortyone consecutive patients were treated with WAL (body-jet®; human med) for mild to moderate body contouring. Patients were given local anesthesia (standardized tumescent solutions) during all three phases of the surgery.

Fat harvesting was accomplished by collecting and separating the aspirated adipose tissue in a sterile container LipoCollector®), without need for washing or centrifugation. Fat grafting by microdroplet technique was performed within two hours of collection. Fat aliquots from five randomly-selected patients were assessed with a trypan blue dye exclusion test within one hour and again six to eight hours after collection."

Results: "A total of 37 females and four males underwent WAL in this series; average body mass index (BMI) was 25.5.

All patients experienced uneventful recovery periods with minimal side effects and no significant complications. Although large volumes of tumescent solution were required during the three phases of the technique, the total volume of infiltration almost equaled the final volume of aspiration. The average infiltration-to aspiration ratio was 1.1 to 1.0 in all cases.

Lidocaine dosage averaged 10.5 mg/kg in Group 1 and 20.0 mg/kg in Group 2. Trypan blue dye exclusion testing indicated that about 90 % of adipocytes expelled the dye after one hour of extraction (resulting in 90 % cell viability). The amount of instilled tumescent fluid, lidocaine dosage, and aspiration volumes appeared to be safe, with minimal blood loss in small and moderate volume liposuction cases."

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H. Peltoniemi, A. Salmi, S. Miettinen et al.: Stem cell enrichment does not warrant a higher graft survival in lipofilling of the breast: A prospective comparative study. Journal of Plastic, Reconstructive & Aesthetic Surgery (2013) 66, 1494-1503.

The prospective comparative study evaluates a "total of 18 women who underwent breast augmentation with water-assisted lipotransfer (WAL)." "Liposuction for retrieval of fat was performed according to the WAL method, as described before, under local anaesthesia and light sedation. The body-jet system was used in combination with the LipoCollector® (Human Med AG, Schwerin, Germany)." "In 10 of the cases, transferred lipoaspirate was enriched with stromal stem cells using the Celution system (Cytori Therapeutics Inc.) Magnetic resonance imaging (MRI)-based volumetric analysis was done preoperatively and 6 months after the procedure.

Results: "MRI volumetry revealed a volume survival of the whole (watery) graft of mean 54 % in the WAL only and of 50 % in the WAL with stem cell-enrichment patients. As centrifugation of the WAL grafts demonstrated an average adipose tissue of 68 %, the average volume survival of adipose tissue itself was 79 % in the WAL only and 74% in the WAL with stem cell-enrichment patients. This difference (4.5 %) was not statistically significant." Conclusion: "We found a high survival rate after WAL and cell enrichment in the presented patients, but not better than in patients purely treated with WAL, without stem cell enrichment. WAL alone is faster (90-150 min less), cheaper (cost of consumables for Celution was over 3.000 euros for each patient), theoretically safer (lower risk of contamination) and offers at least the same take rate. We do not see any advantage in stem cell enrichment by the Celution system."

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This retrospective European multicenter trial included "135 procedures on 28 (35 breasts) postmastectomy patients (mean 52.4 years)." "All women were treated with the water-jet assisted fat grafting method (BEAULITM) combined with additional procedures (NAC reconstruction, contralateral mastoplasty) and evaluated with at least 6 months follow-up (mean 2.6 years). Sonography or mammography, clinical examination, patient questionnaire (10-point Likert scale) and digital photographs were carried out."

Results: "On average the patients received 4 to 6 procedures each with a single volume of 159 ml (±61 ml) over 21 months (range 9 months to 2.5 years). In total 1,020ml (±515 ml) fat were grafted until a complete breast reconstruction was achieved. A complete breast reconstruction with large volume fat grafting is alternatively possible to standard techniques in selected cases. It takes at least 4 to 6 lipotransfers in the course of 2 years. Patients with prior radiotherapy may require even up to 8 sessions over nearly 3 years of treatment." Conclusion: "Our study demonstrates that autologous fat grafting is a suitable way not only to restore local defects after conservative breast reconstructive surgery or radiodamages, but also to achieve complete breast reconstruction."

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T. K. Malan: Breast Augmentation and Reconstruction with Fat Transfer. In: Cosmetic Surgery - Art and Techniques; Shiffman, Melvin A., Di Giuseppe, Alberto (Eds.), 2013, pp 595-603

In this educational book on breast augmentation and reconstruction, the author explains his operative technique as follows: "Following marking the fat is removed from the donor site with the Body-Jet system with Lipocollector to harvest the fat under low pressure (12 in. Hg). Smooth, slow, even passes of harvesting cannula with care taken to not torque the tissue or to rapidly vent cannula causing splash injury as fat rushes into collector. The LipoCollector® drains excess tumescent solution utilizing a 300-micron filter leaving minimally processed fat."

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M. Stabile, K. Ueberreiter, H. E. Schaller, D.L.Hoppe: Jet-assisted fat transfer to the female breast: preliminary experiences. European Journal of Plastic Surgery March 2014

Two specialized European centers. "Retrospective analysis enrolled over 167 female nonsmokers with reconstructive or aesthetic indications treated with at least one procedure according to the standardized protocol of the BEAULITM technique from February 2010 to June 2012. The fat's harvest is carried out in a sterile closed system with a low suction force, avoiding overly mechanical trauma or thermal damage. The WAL system uses nearly 70 % less tumescent fluid than other liposuction systems. In this way, local swelling and bloating are minimized, preserving the contours of the target area. With the LipoCollectorTM system, neither centrifugation nor additional washing are required."

Results: "The study included. 487 autologous jetassisted fat transfer procedures with minor complications (5.35 %) like tiny oil cysts formations or hematoma of the donor site.

Low postoperative pain has been reported during the first week, and final aesthetic evaluation showed good to excellent results. We observed a higher frequency of procedures, especially in irradiated patients (>3.62) compared to the rest of reconstructive cases (>2.78)."

Conclusions: "Water-jet-assisted liposuction with consecutive immediate mammary fat injection is a procedure with a short hospitalization and low complication rate."