# hydro surgery

## news

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# „De novo“ breast reconstruction with autologous fat grafting in postmastectomy patients\*

D.L. Hoppe<sup>1</sup>, K. Ueberreiter<sup>1</sup>, S. Kauhanen<sup>2</sup>, Y. Surlemont<sup>3</sup>

<sup>1</sup>Park - Klinik Birkenwerder, Birkenwerder, Germany, <sup>2</sup>Helsinki University Hospital, Plastic and Reconstructive Surgery, Helsinki, Finland, <sup>3</sup>Clinique Saint Antoine, Chirurgie plastique et esthétique, Rouen, France

## Introduction

„Fat tissue transplant is been approved as a safe and reliable method for breast reconstructive issues in oncologic patients. Various clinical applications for autologous fat grafting to the breast have been described in literature (1). Correction of soft tissue defects and deformities after breast cancer or insufficient reconstruction can be achieved (2). However complete breast reconstruction after total mastectomy is mentioned in few cases (3,4). Our multicenter study investigates in larger series the aspect of lipotransfer for total breast reconstruction after primary ablative surgery.

## Material/Methods

In a retrospective clinical trial we analysed the data of 3 centers in Europe including a total of 54 patients (62 breasts) after uni- or bilateral radical mastectomy. Over 300 procedures were performed between 2008 and 2012. Inclusion criteria were fulfilled by 21 patients. Autologous fat grafting was done according to standardised protocol of the BEAULI Method (4). Clinical outcomes, patient satisfaction and aesthetic results were evaluated after a 6 months follow-up period. If available, comparative MRI volumetry was done for quantification (5).

## Results

All patients showed remarkable increase of subcutaneous fat. On average 3,8 procedures with a single median volume of 150 ml ( $\pm$  57,4 ml) were performed till the end of treatment. A lipotransfer of 600 ml ( $\pm$  177 ml) was required for complete breast reconstruction. The median operation time was 40 min, but ranged to 60 min ( $\pm$  29 min) when additional operative steps

(NAC reconstruction, abdominal advancement, contralateral mammoplasty) were carried out. 44,7 % of the patients underwent their intervention in analgesic sedation. As a postoperative complication small singular oilcysts were estimate in 4 (2,6 %) cases. Once granuloma was detached, but no infections. The majority of patients was very satisfied with the final result.

## Discussion

Breast lipografting in oncological patients is still considered as a delicate and risky surgical approach. However this study supports the use of standardized protocols for autologous lipoaspirate grafting with longterm experience in order to provide optimized cosmetic outcome. It emphasizes, that a „de novo“ breast reconstruction is possible within 3-5 procedures of BEAULI Method. These are combinable simultaneously with other oncoplasty techniques. Our reconstructive treatment is gradually attained after 12 -18 months at least. Patients, who declined immediate reconstructive, prothesis - based approaches or experienced numerous complications after breast surgery are very satisfied with BEAULI Method. In 4 years of experience we notice an increasing demand. Although further studies are certainly required to determine definitive guidelines for lipotransfer in breast surgery, it might be regarded as an evolving and promising alternative to conservative methods.“


## References

- (1) Coleman SR, Saboeiro AP. Fat grafting to the breast revisited: safety and efficacy. *Plast Reconstr Surg.* 2007;119 (3): 775-85.
- (2) Rigotti G et al. Determining the Oncological Risk of Autologous Lipoaspirate

Grafting for Post-Mastectomy Breast Reconstruction. *Aesth Plast Surg* 2010; 34: 475-480.

- (3) Babovic S. Complete breast reconstruction with autologous fat graft - a case report. *J Plast Reconstr Aesthet Surg.* 2010; 63 (7): e561-3.
- (4) Illouz Y.G., Sterodimas A. Autologous Fat Transplantation to the Breast: A Personal Technique with 25 Years of Experience. *Aesth Plast Surg* 2009; 33:706-715
- (5) Ueberreiter K et al. BEAULI - A New and Easy Method for Large Volume Fat Grafts. *Handchir Mikrochir Plast Chir* 2010; 42: 379 - 385.
- (6) Herold C, Ueberreiter K et al. Is There a Need for Intrapectoral Injection in Autologous Fat Transplantation to the Breast? - An MRI Volumetric Study. *Handchir Mikrochir Plast Chir.* 2011; 43: 119 - 124.

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## Does stem cell enrichment enhance fat uptake in jet-assisted fat transfer for breast augmentation?\*

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„Body-Jet-assisted fat transfer for augmentation and reconstruction of breasts is popular in Finland. Overall, more than 100 breast augmentation procedures by lipotransfer (BEAULI method) have been performed by the author with mean follow-up of 18,9 months (4-36). The purpose of the study was to determine if enriching the transferred fat with autologous stem cells will enhance uptake of BEAULI fat. In a prospective study, eleven women had breast augmentation with enrichment with Cytori Celution method and seven women had non-enriched, normal BEAULI transfer.

In volumetric analysis based on MRI, no superiority of either method could be determined. Individual differences in especially recipient area (elasticity, thickness of subcutaneous fat layer) seem to be of more importance than the amount of stem cells in transferred fat. Healthy breasts contain billions of stem cells, and atraumatically collected fat contains millions, so enrichment may not be necessary or even cost beneficial in BEAULI fat transfer. In scarred or damaged stem cell poor tissue, enriching might be beneficial.“

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„Unprepared samples are more viable than the samples that were washed or centrifuged...The less one manipulates the fat graft and the more quickly it is reinjected, the higher the chances of survival are.“ P. Smith et al.: Autologous Human Fat Grafting: Effect of Harvesting and Preparation Techniques in Adipocyte Graft Survival. *PRS*, 117: 1836-44 (2006)

„Presence of greater amounts of pre-adipocytes in the non-centrifuged adipose tissue cultures and more distinctly expressed cell proliferation. - Our data suggests that with transplantation of non-centrifuged adipose tissue more active pre-adipocytes are applied which could possibly lead to better potential chances of survival and even de novo development of fat.“ Khater et al.: Clinical and Experimental Study of Autologous Fat Grafting After Processing by Centrifugation and Serum Lavage. *Aesth Plast Surg* 33:37-43 (2009).

Cell counts per high-powered fields of intact adipocytes and nucleated adipocytes and adipocyte cross-sectional area were significantly greater in samples processed by sedimentation, compared

with those by centrifuging and washing“.

Rose et al.: Histologic Comparison of Autologous Fat Processing Methods. *Opht Plastic & Reconstr Surg.* Vol. 22, 3:195-200 (2006).

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Von Professor M.C. Harmsen University of Groningen, Cardiovascular Regenerative Medicine Research Group, Department of Pathology & Medical Biology, University Medical Centre Groningen, Groningen, Netherlands

„The adventitia of the vasculature of adipose tissue harbours mesenchymal stem cells with potent regenerative capacity, the so-called Adipose Tissue-derived Stem Cells (ADSC). In practice, ADSC are isolated as the plastic-adherent cells in the stromal vascular fraction of fat tissue. Their frequency is around 100x10<sup>6</sup> per litre of lipoaspirated fat. ADSC can differentiate into multiple lineages including fat, bone, cartilage, smooth muscle cells and pericytes. Yet, their main virtue resides in the secretion of a plethora of paracrine factors such as cytokines, chemokines,

growth factors and non-proteinaceous factors that greatly augment healing of damaged tissues. ADSC have been successfully employed in various regenerative medicine-based strategies that range from tissue-engineering of bone and cartilage to cardiovascular repair. Plastic surgeons have long noted that lipotransfers reduce local scars and rejuvenate skin. We hypothesized that the ADSC present in fat, have scar-reducing i.e.anti-fibrotic potential. We tested our hypothesis by subjecting differentiating (i.e.TGFβ treated) human dermal fibroblasts (HDF) to conditio-

ned culture media of ADSC (ADSCcm). We showed that ADSCcm not only suppressed the TGFβ-driven differentiation of HDF, but that also their contractile activity was reduced significantly. Furthermore, ADSCcm reduced the production and secretion of extracellular matrix by HDF, while the secretion of active matrix remodelling enzymes (MMP1, 2, 3 and 9) was upregulated. We conclude that ADSC may be suitable therapeutic cells to harness (hypertrophic) scarring.“

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(\*Patient: female, 27 years; liposuction from thighs; body-jet®evo RANGE 2 during infiltration, RANGE 1 during suction; vacuum -500 mbar) Source: University of Rostock, Department of Cell Biology; results by FACS analysis: Seracell Lab, Instrument SN: AM07032 SW VersionCXP v2.2

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(\*Patient: female, 59 years; liposuction from abdomen; body-jet®evo RANGE 2 during infiltration, RANGE 1 during suction; suction vacuum -500 mbar; WAL cannula 3.8 rapid) Source: University of Rostock, Department of Cell Biology; results by FACS analysis: Seracell Lab, Instrument SN: AM07032 SW VersionCXP v2.2

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### Ein neuer Fokus der human med AG